

SALL4 (6E3)

Format	Catalog No.	Pack size	Dilution
Concentrated	G384 A, B, C	0.1, 0.5, 1.0 mL	1:100
Prediluted	G384 AA, BB	6.0, 3.0 mL	Ready to use

SPECIES: Mouse monoclonal

IMMUNOGEN: SALL4 (954-1054)

CLONE: 6E3

ISOTYPE: IgG1/ kappa

SUPPLIED AS – Buffer with protein carrier and preservative

BACKGROUND –

Through its modulation of Oct4, SALL4 is necessary for the preservation of embryonic stem cell pluripotency. SALL4 has been shown to be a sensitive and specific marker for ovarian primitive germ-cell tumours and seminomas. Research has shown that high SALL4 staining is present in more than 90% of tumour cells in embryonal carcinomas and intratubular germ-cell neoplasias. Furthermore, all 31 yolk sac tumours (5 paediatric and 26 postpubertal) had tumour cells that were strongly positive for SALL4 labelling but negative for Oct4. Research indicates that SALL4 antibody is a more promising pan germ-cell marker than Oct4 and PLAP antibodies.

SPECIES REACTIVITY : Humans; others not tested.

POSITIVE CONTROL : Seminoma

CELLULAR LOCALIZATION : Nuclear

TITER/WORKING DILUTION : Immunohistochemistry (Frozen and Formalin-fixed): 1-2ug/ml

Flow Cytometry: 1-2ug/million cells

Immunofluorescence: 1-3ug/ml

Western Blotting: 2-4ug/ml

MICROBIOLOGICAL STATE : This product is not sterile.

LIMITATIONS AND USES:

1. For Research Purpose Only (RUO)
2. Not to be consumed internally.
3. This product is designed for qualitative immunohistochemistry using formalin-fixed, paraffin-embedded
4. tissue sections from both normal and malignant tissues, which will be examined under a light microscope.
5. If the reagent gets hazy, don't use it.
6. Never use after the expiration date.
7. When working with reagents, exercise caution.
8. Not sterile.



4448 Ammendale Rd Beltsville, MD
20705, USA



www.genebiosolution.com



info@genebiosolution.com



+1 (408) 580-1396

MATERIALS REQUIRED BUT NOT PROVIDED

1. Positive Tissue Control: Routinely processed, neutral buffered formalin-fixed, paraffin-embedded Bladder or urothelial carcinoma
2. Negative tissue control (internal or external)
3. Microscope slides and coverslips
4. Staining jars or baths
5. Timer Instructions for Use (IFU)
6. Xylene or xylene substitute
7. Ethanol or reagent alcohol
8. Deionized or distilled water
9. Heating equipment or enzyme for tissue pretreatment step
10. Detection system
11. Chromogen
12. Wash Buffer
13. Hematoxylin
14. Antibody diluents
15. Peroxide Block
16. Light Microscope
17. Mounting medium
18. Avidin-Biotin Blocking Reagents for use with streptavidin biotin detection

PROCEDURE –

1. **Tissue Section Pretreatment (Necessary):** Citrate Plus pretreatment greatly improves the staining of tissue slices that have been paraffin embedded and fixed with formalin.
2. **Incubation Time for Primary Antibody:** We recommend 30 minutes at room temperature for incubation.
However, the user should choose the best incubation time based on the staining system used and the fixing circumstances.
3. **Visualisation:** To get the highest staining intensity, we advise combining the "DAB Chromogen/Substrate Bulk Pack (High Contrast)" with the "HRP Pack".

PRECAUTIONS –

1. Contains 0.09% w/v sodium azide as a preservative.
2. Avoid pipetting by mouth.
3. Reagents and specimens should not come into touch with skin or mucous membranes.
4. Prevent reagent contamination by microorganisms, as this could lead to an increase in nonspecific staining.
5. There are no dangerous materials in this product.

TECHNICAL SUPPORT

For technical assistance, please contact Genebio Solution's Technical Support at www.genebiosolution.com



4448 Ammendale Rd Beltsville, MD
20705, USA



www.genebiosolution.com



info@genebiosolution.com



+1 (408) 580-1396