

GBS-SP001 HIGH PRECISION LASER GLASS SLIDE PRINTER



Products Features:

01. Super large capacity design. The number of slices for a single baking is no less than 220. Set floating, baking, baking in one, hot air circulation design;
02. The spreading box has the function of underwater light transmission, and the spreading sheet is not affected by external light, providing the best position to observe the specimen and the spreading sheet;
03. Dual heating, sensing, and control systems. When one set of systems malfunctions, it automatically switches to the other set. The normal operation of the instrument is not affected;
04. Large screen LCD display;
05. With the function of timing startup and shutdown, you can set the timing startup and shutdown time, and then set the day of the week to start and the day of the week not to start, the stall baking machine can be no longer manual management, automatic switch on and off every day;
06. PID fuzzy control technology, effectively overcome thermal inertia, constant and accurate temperature control;
07. Has a separate baking sheet, baking sheet opening key, when only the spreading function can choose to open or close baking sheet, baking sheet, baking sheet double timing;
08. Baking sheet fast, uniform, good effect;
09. Fully reflects the ergonomic Angle of 60 degrees baking sheet insert, can bake ≥ 100 pieces;
10. Two level safety protection, original heating components;
11. Spreading temperature: room temperature $\sim 99^\circ\text{C}$ adjustable;
12. Baking temperature: room temperature $\sim 99^\circ\text{C}$ adjustable;

Technical Parameter

Item	Parameter
Display Screen	LCD Display
Maximum Quantity per Drying	220 PCS
Spreading Temperature	Room Temperature $\sim 99^\circ\text{C}$ adjustable
Baking Temperature	Room Temperature $\sim 99^\circ\text{C}$ adjustable
Voltage / Frequency	220V / 50HZ
Power	600W
Dimension	330MM x 430MM x 150MM
Weight	7KGS

